

# Mechanical Comparison of UP and PETG Polymers

---

**CM000058880**

Andrea Westlie, Unplastic

**Authors**      Stephanie Currie, M.S.  
Josh Campbell, Ph.D.

**Reviewers**    Shivesh Sivakumar, M.S.


**April 15, 2026**

## Goal:

- Evaluate mechanical behavior of PETG and UP polymers under conditions simulating clinical use.

## Methods:

| Analysis                 | Instruments & Operating Parameters                                                                                                                                                                                                                                          |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Flex Testing (ASTM D790) | <div>Instron 68SC-5 Universal Testing Machine (<i>Appendix 2</i>)</div> <ul style="list-style-type: none"><li>3-point bend fixture w/ 5-mm radius anvils</li><li>ASTM D790 Procedure A, Type I*</li><li>500 N load cell</li><li>Ambient Environment: 22°C, 47% RH</li></ul> |
| Creep Testing            | <div>Anton Paar MCR 702 Multidrive (<i>Appendix 3</i>)</div> <ul style="list-style-type: none"><li>Constant tensile load of 1.0 N</li><li>Extensional strain monitored for 60 minutes</li><li>Ambient temperature and 37°C</li></ul>                                        |



Flexural Samples – After Testing

## Summary of Results:

- Flexural testing shows PETG to have higher Yield Strength and Flexural Strength, while UP has higher Modulus of Elasticity.
- Creep testing shows UP polymer to be more resistant to creep deformation at both ambient and 37°C.
- Stress relaxation test was not feasible using Anton Paar MCR 702.
  - An alternative test method is proposed in “Next Steps”.

\* Strain rate of 0.01 mm/mm/min (Procedure A) with strain calculated from crosshead displacement (Type I)



| Specimen         | Modulus of Elasticity* [GPa] | Yield Stress (0.2% offset) [MPa] | Flexural Strength [MPa]** |
|------------------|------------------------------|----------------------------------|---------------------------|
| PETG - 1         | 2.73                         | 92.0                             | 96.9                      |
| PETG - 2         | 2.73                         | 87.5                             | 96.0                      |
| PETG - 3         | 2.70                         | 90.5                             | 95.7                      |
| <b>Average</b>   | 2.72                         | 90.0                             | 96.2                      |
| <b>Std. Dev.</b> | 0.02                         | 2.3                              | 0.1                       |

| Specimen         | Modulus of Elasticity* [GPa] | Yield Stress (0.2% offset) [MPa] | Flexural Strength [MPa] |
|------------------|------------------------------|----------------------------------|-------------------------|
| UP- 1            | 3.13                         | 51.8                             | 77.6                    |
| UP - 2           | 3.02                         | 43.9                             | 71.6                    |
| UP - 3           | 3.05                         | 74.1                             | 82.6                    |
| <b>Average</b>   | 3.07                         | 56.6                             | 77.3                    |
| <b>Std. Dev.</b> | 0.06                         | 15.7                             | 5.5                     |



\* 0.5 – 1.0% Chord Modulus, per ASTM D790, section 5.1.8.3.

\*\* Per ASTM D790 requirements, test terminates at 5% strain.

- PETG samples show higher Yield stress and Flexural Strength.
- UP shows higher Modulus of Elasticity.
- Deformation behavior is more variable with the UP samples.

- As indicated by higher standard deviations for Yield Stress & Flexural Strength.

# Creep Testing Results



- Extensional Strain as a function of time (left plot) shows that UP polymer has higher strains at the applied load of 1.0N.
- Looking at the change in strain from the start of the test (right plot) shows that UP is more dimensionally stable with time at both room temperature and 37 °C.
- Both polymers show increased creep strain as temperature is increased (for the temp. range investigated).

- Stress Relaxation Testing was attempted on the Anton Paar MCR 702 (Appendix 3), by applying a fixed strain and monitoring stress over a 24-hour period.
  - This instrument was unable to sufficiently hold the prescribed strain value.
- The Instron 68SC-5 universal testing machine (Appendix 2) will likely be better able to hold the prescribed strain value.
  - Covalent proposed attempting the stress relaxation test on this instrument.
  - Due to lack of an environmental chamber, it will only be possible to perform this test at ambient temperature.

## Flexural Testing

| Specimen | Width (mm) | Thickness (mm) | Load Span (mm) |
|----------|------------|----------------|----------------|
| PETG 1   | 12.83      | 1.37           | 22             |
| PETG 2   | 12.88      | 1.37           | 22             |
| PETG 3   | 12.86      | 1.37           | 22             |
| UP 1     | 12.94      | 0.86           | 14             |
| UP 2     | 12.91      | 0.84           | 14             |
| UP 3     | 12.94      | 0.89           | 14             |

## Creep Testing

| Specimen | Width (mm) | Thickness (mm) | Length (mm) |
|----------|------------|----------------|-------------|
| PETG, RT | 9.99       | 1.37           | 24.95       |
| PETG, 37 | 9.79       | 1.37           | 19.08       |
| UP, RT   | 10.17      | 0.76           | 15.90       |
| UP, 37   | 9.40       | 0.68           | 15.79       |




| Parameter        | Value                                                  |
|------------------|--------------------------------------------------------|
| Maximum Load     | 5000 N                                                 |
| Load Accuracy    | ± 0.5% of reading down to 1/1000 of load cell capacity |
| Speed Range      | 0.005 to 2000 mm/min                                   |
| Crosshead Travel | Up to 1120 mm (44 in.)                                 |
| Software         | Bluehill Universal                                     |

# Appendix 3 – Instrument: Anton Paar MCR 702 Rheometer

DMA is a measurement of a material’s viscoelastic properties as a function of time, temperature, and frequency.

E’, E”, tan delta, frequency sweep, temperature sweep, time-temperature superposition (TTS).



| Mode            | Length (mm) | Width (mm) | Thickness (mm)                | Temperature range (°C) | Stiffness range (MPa) | Example materials                                  |
|-----------------|-------------|------------|-------------------------------|------------------------|-----------------------|----------------------------------------------------|
| Tension/Torsion | 5-35        | 5-15       | 0.01 - 5                      | -100 - 350             | 0.1 - 10000           | Films, fibers, flat samples                        |
| Compression     | 1-40        | 1-40       | 0.5 - 5                       | -100 - 220             | 0.0001 - 100          | Foams, elastomers, cylinders/discs                 |
| 3-pt Bend       | 45-50       | 10-12      | stiff = 0.2-0.3<br>soft = 1-2 | -100 - 350             | 100 - 1000000         | Composites, glass, ceramics, cured resins/polymers |