

Built to Handle Everyday Exposure

Tensile performance after repeated cleaning and oral-bacteria biofilm exposure

Study 1 | Tensile Testing After Cleaning - Covalent Material Insights, May 2026

Study 2 | *S. mutans* Biofilm Degradation Pilot - University at Buffalo Materials Testing Facility, Apr/May 2026

Staining | 24- and 72-hour immersion in common beverages

Unplastic | Summary followed by both full test reports

Summary of Results

Both full reports follow

STUDY 1 - COVALENT MATERIAL INSIGHTS

Repeated cleaning

ASTM D638 Type V dogbones, Instron 68SC-5 with AVE3 video extensometer, n = 3 per group.

Group	Modulus [GPa]	UTS [MPa]	Strain at UTS [%]
CTRL	2.48 ± 0.13	37.02 ± 0.49	4.54 ± 0.46
CS14d	2.82 ± 0.53	37.10 ± 0.82	4.49 ± 0.29
UP28	2.66 ± 0.24	36.55 ± 0.70	4.85 ± 0.20

Result: All three groups showed comparable tensile behavior, with mean ± SD ranges overlapping for every property. Tensile strength stayed within 36.5-37.1 MPa across groups.

Values are mean ± standard deviation as reported by Covalent (report CM000060037, May 12, 2026).

STUDY 2 - UNIVERSITY AT BUFFALO (PILOT)

Oral-bacteria biofilm

Moses PHA, MosesPHA-2, and a PETG aligner control aged 2 weeks in *S. mutans* NCTC 10449 (BHI + 1% sucrose), deionized water, or air; cleaned with 70% IPA between cycles. n = 2 per condition.


Max tensile strength [MPa]	Air	Water	Biofilm
PETG	56.54	55.40	55.23
Moses PHA (PB3410)	41.19	26.08	36.21
MosesPHA-2 (PV3020)	37.48	35.74	35.88

Lab conclusions: No noticeable aesthetic changes to any specimen; biofilm grew mainly on roughened surfaces; the pilot showed no discernible difference in strength after biofouling exposure.

Values as reported by the lab. Pilot study with n = 2, so no statistics were reported. In the full report, charts label Moses PHA as PB3410 and MosesPHA-2 as PV3020.

Stain Resistance

24- and 72-hour beverage immersion



IMMERSION TEST

Dogbone specimens were immersed in six everyday beverages and photographed after **24 hours** (left column) and **72 hours** (right column).

Beverages

Red wine
Orange juice
Diet Coke
Celsius
Black instant coffee
Black tea

The vertical specimen at far left is shown for comparison.

Visual comparison only; photo shown unedited.

STUDY 1 - FULL REPORT

Tensile Testing After Cleaning

Covalent Material Insights | ASTM D638, Type V | May 2026

Tensile Testing of Polymer Samples per ASTM D638, Type V

CM000060037

Andrea Westlie, Unplastic

Authors Josh Campbell, Ph.D.

Reviewers Shivesh Sivakumar, M.S.

May 12, 2026

Goal:

- Quantify polymer performance via tensile testing per ASTM D638.

Methods:

- Mechanical Testing with Instron 68SC-5.
- Strain measurement via non-contact Instron AVE3 video extensometer.
- Method details in Appendix.


Summary of Results:

- The 3 polymer materials showed very comparable material properties.



Tensile Properties (Average $\pm$ Std. Dev.*)

Material	Modulus of Elasticity [GPa]	Tensile Strength, UTS [MPa]	Strain at UTS [%]
CTRL	2.48 $\pm$ 0.13	37.02 $\pm$ 0.49	4.54 $\pm$ 0.46
CS14d	2.82 $\pm$ 0.53	37.10 $\pm$ 0.82	4.49 $\pm$ 0.29
UP28	2.66 $\pm$ 0.24	36.55 $\pm$ 0.70	4.85 $\pm$ 0.20

* $n=3$



Results – All Material Samples




Results Table (previously presented on slide 2)

Material	Modulus of Elasticity [GPa]	Tensile Strength, UTS [MPa]	Strain at UTS [%]
CTRL	2.48 ± 0.13	37.02 ± 0.49	4.54 ± 0.46
CS14d	2.82 ± 0.53	37.10 ± 0.82	4.49 ± 0.29
UP28	2.66 ± 0.24	36.55 ± 0.70	4.85 ± 0.20

- All materials showed comparable tensile behavior with Mean +/- Standard Deviation range overlapping for all properties.
- Tensile data are plotted up to the Ultimate Tensile Strength.
 - Strain after UTS was minimal.
 - Artifacts in strain signal occurred at fracture due to elastic rebound being high relative to the total strain values.
- Due to above factors, material ductility is compared in terms of Strain at the Ultimate Tensile Strength.


Results – CTRL Sample



Replicate	Modulus of Elasticity [GPa]	Tensile Strength, UTS [MPa]	Strain at UTS [%]	Thickness [mm]
1	2.47	37.44	5.07	0.78
2	2.36	36.49	4.22	0.85
3	2.62	37.13	4.33	0.86
Mean	2.48	37.02	4.54	0.83
Std. Dev.	0.13	0.49	0.46	0.04

- Highly repeatable tensile behavior for 3 replicates.


Results – CS14d Sample



Replicate	Modulus of Elasticity [GPa]	Tensile Strength, UTS [MPa]	Strain at UTS [%]	Thickness [mm]
1	2.65	37.86	4.75	0.88
2	2.40	37.19	4.56	0.82
3	3.42	36.24	4.18	0.78
Mean	2.82	37.10	4.49	0.83
Std. Dev.	0.53	0.82	0.29	.05

Slight differences, as compared to CRTL specimen:

- Higher Modulus mean and standard deviation.
- Higher variation in Ultimate Tensile Strength (UTS).
- Lower variation in Strain at UTS.



Replicate	Modulus of Elasticity [GPa]	Tensile Strength, UTS [MPa]	Strain at UTS [%]	Thickness [mm]
1	2.76	35.90	4.69	0.74
2	2.39	36.45	5.08	0.87
3	2.84	37.30	4.79	0.65
Mean	2.66	36.55	4.85	0.75
Std. Dev.	0.24	0.70	0.20	0.11

- Slight differences, as compared to CRTL specimen:
- Higher Modulus and Strain at Ultimate Tensile Strength (UTS).
 - Lower Ultimate Tensile Strength.

Instrument	Instron 68SC-5 (Appendix 2)
Test Methods	• ASTM D638
Samples	• Type V dog bone samples, received ‘ready to test’
Load Train	• 500 N load cell • 250 N screw action grips • 25 mm x 25 mm serrated jaw faces
Testing Rate	• 1 mm/min
Testing environment	• Ambient: 24.5 °C, 41.8% RH
Extensometer	• Instron AVE3 (<i>Appendix 3</i>)



Parameter	Value
Maximum Load	5000 N
Load Accuracy	± 0.5% of reading down to 1/1000 of load cell capacity
Speed Range	0.005 to 2000 mm/min
Crosshead Travel	Up to 1120 mm (44 in.)
Software	Bluehill Universal



Parameter	Value
Applications	Plastics, metals, composites, textiles, films, elastomers, paper, components, biomaterials (nearly any material)
Loading	Tensile, compression, bend
Resolution	0.5 μm
Accuracy	$\pm 1 \mu\text{m}$ or 0.5% of reading (whichever is greater)
Gage length range	4 – 175 mm

STUDY 2 - FULL REPORT

S. mutans Biofilm Degradation Pilot

University at Buffalo Materials Testing Facility | April/May 2026

Unplastic Tensile S.Mutans Biofilm Degradation Pilot Results

April/May 2026

Materials Testing Facility

University at Buffalo

Project Scope


- Artificial aging of two experimental (PHA) and one PETG aligner groups dogbone specimens in S.Mutans (2ea), Deionized water (2ea), and air (2ea) for two weeks
- Tensile testing to failure of specimens after two weeks

Protocol Description


- Based on ASTM D638-22
 - Type V dogbone specimens
- 70% isopropanol wipes cleaning to disinfect and remove biofilms from specimens before each biofilm cycle (5 cycles; every 2-3 days; 3 weeks)
 - S.Mutans NCTC 10449 strain in brain heart infusion broth with 1% sucrose
- Tensile testing
 - Specimen width and thickness measured at three locations along gage length to calculate respective average dimensions
 - Wedge grips used to fixture specimens to universal testing machine (Dillion Quantrol TC2i)
 - Crosshead speed: 10mm/min until break
 - Elastic slope recorded (Young's Modulus calculated)
 - Nominal Displacement recorded (nominal strain calculated)
 - Maximum tensile load recorded (max tensile strength calculated for failure at yield or fracture)

Tensile Testing Results

Young's Modulus




Maximum Tensile Strength at Yield





Representative Curves and Specimens




PETG - Air




PETG - Water




PETG – S.Mutans





Moses PHA - Air




MosesPHA - Water




MosesPHA - S.Mutans






MosesPHA-2 - Air



MosesPHA-2 - Water



MosesPHA-2 – S.Mutans



Conclusions and Recommendations

- No noticeable aesthetic changes to any specimens
- Biofilms appeared to mainly grow on roughened surfaces
- The pilot showed no discernable difference in strength of the specimens after exposure to biofouling